# HD 8602

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ENGLISH

## HD 8602

- 
- 1 8-pole DIN connector for the TP 870 temperature probe input.
- 
- 2 Low battery indication.
- 
- 3 HOLD symbol.
- 
- 4 — symbol.
- 
- 5  $\Delta$  symbol.
- 
- 6  key. Blocks the reading while internally the instrument continues updating data.
- 
- 7  key. Sets pH measurement or if the instrument is in the calibration mode this activates calibration of the asymmetry potential.
- 
- 8  key. Increases the value shown on the display during the parameters setting.
- 
- 9  key. When this key is pressed the instrument is set for carrying out or confirming calibration.
- 
- 10 BNC socket for connecting the electrode.
- 
- 11 Display.
- 
- 12  $^{\circ}\text{C}$  and  $^{\circ}\text{F}$  symbols.
- 
- 13 mV symbol.
- 
- 14 pH symbol.
- 
- 15  key. Switches the instrument On or Off.
- 
- 16  key. Selects reading in  $^{\circ}\text{C}$  or  $^{\circ}\text{F}$ .
- 
- 17  key. Selects mV measurement; if the instrument is in the calibration mode this activates slope calibration.
- 
- 18  key. Decreases the value shown on the display during the parameters setting.
-

## 1. FOREWORD

The pH concept was introduced in 1909 by S.P.L. Sorenson. It is defined as:

$$(1.1) \quad \text{pH} = -\log_{10} C_{\text{H}^+}$$

where  $C_{\text{H}^+}$  is the concentration (based on weight/volume) of hydrogen ions  $\text{H}^+$ . The pH is therefore a unit for measuring the concentration of  $\text{H}^+$  ions in a watery liquid, just as the metre is a unit for measuring length and the kilo for mass.

G.N. Lewis discovered that ionic activity does not coincide with concentration, consequently a pH electrode measures ionic activity and not the concentration of the  $\text{H}^+$  ion. Based on the connections between ionic activity and concentration, the pH value measured with a pH electrode deviates from the real value of the concentration of  $\text{H}^+$  ions, the higher the ionic concentration of the solution to be analysed.

At  $25^\circ\text{C}$  pure water has a pH value of 7 (neutral). Solutions with a pH value  $< 7$  are called acid solutions; solutions with  $\text{pH} > 7$  are called alkaline solutions. Table 1.1 shows the pH values of some substances in everyday use.

**Table 1.1**

pH of various substances		
	0	
	1	← Sulphuric acid 4.9% (0.1 N)
	2	← Hydrochloric acid 0.37% (0.1 N)
Lemon juice →	3	
Orange juice →	4	← Acetic acid 0.6% (0.1 N)
Beer →	5	
Cheese →	6	← Prussic acid 0.27% (0.1 N)
Milk →	7	
Pure water →	8	
Egg white →	9	← Bicarbonate of soda 0.84% (0.1 N)
Borax →	10	← Potassium acetate 0.98% (0.1 N)
Milk of magnesia →	11	← Ammonia 0.017% (0.01 N)
	12	← Ammonia 1.7% (1 N)
	13	← Caustic soda 0.04% (0.01 N)
	14	← Caustic soda 4% (1.0 N)

Some acids dissociate strongly in water, forming ions with a positive charge and ions with a negative charge. An acid with a high degree of dissociation, for example sulphuric acid, is a strong acid. Most chemical reactions in aqueous solutions are strongly influenced by the pH value, as regards both speed and final equilibrium. The  $\text{H}^+$  ion may act only as a catalyst, or it may take an active part in the reaction.

## 2. PRINCIPLE OF MEASUREMENT

At present, two methods are used for measuring pH. The first is a colorimetric method based on substances which change colour in certain acid or alkaline environments. This category includes pH revealing papers which change colour according to the pH of the solution. The degree of precision offered is fairly low; a precision of  $\pm 0.1$  pH may be obtained only in a laboratory on account of the intrinsic difficulties of this method:

- a) the addition of the revealing agent may alter the pH value;
- b) the change in colour may be concealed by the presence of suspended particles;
- c) precision depends on the colour of the solution.

The second technique for measuring pH is the potentiometric method, which has been universally adopted by industry. pH is assessed by measuring the difference of potential (voltage) between two electrodes immersed in the solution that is to be analysed. The two electrodes may be physically separate, in which case there is a galvanic cell formed by the glass electrode, the solution with unknown pH and the reference electrode. For ease of handling, the glass electrode and the reference electrode are frequently enclosed in a single tube, in which case it is called a combined electrode. The glass electrode is composed of a glass tube (fig. 2.1) at the bottom end of which is welded a special pH-sensitive glass membrane 9. The inside of the tube is filled with a buffer solution 7 in which the electrode 8 is immersed.

To establish the galvanic potential, a high concentration of KCl is added, for example 3.0 moles, 3.5 moles or a saturated solution of KCl. At the top end of the glass tube is the head 5 through which passes a cable 4 with screen 3, insulation 2 and core 1 (in contact with the electrode 7). In place of the head shown in figure 2.1 there may also be a connector. All the connections (connectors with screened cables) must have a high degree of insulation on account of the high internal impedance of the glass electrode (30...1000M $\Omega$ ).

The reference electrode may be compared to the cold junction when measuring temperature by thermocouple. By measuring the difference of potential between the glass electrode and the reference electrode, we obtain a voltage proportional to pH. In order to obtain reproducible values, the potential generated by the reference electrode must be stable through time and must not present any hysteresis phenomena in the event of temperature variations. Two main types of reference electrode are used today:

- a) the calomel electrode which has the advantage of excellent reproducibility and stability, but it can be used only up to about 70°C.
- b) the silver or silver-chloride electrode can be used with temperatures of over 100°C. However, the reproducibility of the semi-cell potential may present a fair degree of hysteresis in the event of temperature variations.

Physically, an Ag/AgCl reference electrode is composed of a glass or plastic tube in which is applied a porous septum made of a ceramic disc or plaited platinum wires. The inside of the tube is filled with a solution or gel of AgCl and KCl in which is immersed a silver wire covered with AgCl. The function of the porous septum is to create a liquid joint with the solution that is to be analysed.

A small but constant flow of electrolyte from the inside of the electrode to the outside prevents polluting substances from spreading into the reference solution, which would lead to variations of the semi-cell potential.

The necessity of a continuous flow creates the following problems:

## GLASS ELECTRODE


Fig. 2.1

- a) The porous septum may become blocked on account of crystallization of the electrolyte inside. This phenomenon is particularly pronounced at low temperatures.
- b) If the solution to be analysed is under pressure, a slightly higher pressure must be created inside the electrode.
- c) There is a continuous loss of electrolyte.

The glass electrode follows Nernst's law fairly closely, according to which for a variation of pH there is a corresponding signal variation of 59.16 mV (at 25°C). Nernst's voltage also depends linearly on temperature, therefore the instrument must be equipped with some form of manual and/or automatic compensation (as provided in the HD 8602 pH-meter).

### 3. BUFFER SOLUTIONS

The precision of the pH measurements depends on the precision with which the pH of the buffer solutions is known. The National Bureau of Standards (NBS) has measured the pH of a number of buffer solutions and established their value up to 3 decimal places. One must, however, bear in mind that the third decimal place indicates only the reproducibility of the measurements and not absolute precision, which is estimated as  $\pm 0.01$  pH.

Obviously the standard solutions used for calibrating pH-metres must be very stable, consequently we only use solutions whose pH value is not much influenced by the addition of acids and bases or of water (that is, by dilution). These properties are to be found in solutions containing a weak acid and a salt of the same, or a weak base and a salt of the same.


**Table 1 - Buffer Solutions**

Name	Ingredients in 1000 cm <sup>3</sup> H <sub>2</sub> O	pH at 25°C(*)	Degree of buffering $\beta$	Effect of dilution $\Delta$ pH
Potassium quadroxalate	12.61 g KH <sub>3</sub> (C <sub>2</sub> O <sub>4</sub> ) <sub>2</sub> • 2 H <sub>2</sub> O	1.70	0.11	
Acid potassium phthalate	10.12 g KH C <sub>8</sub> H <sub>4</sub> O <sub>4</sub> , dried for 2 hours at 110...130°C	4.01	0.016	+0.04
Phosphate	3.63 g KH <sub>2</sub> PO <sub>4</sub> + 5.68 g Na <sub>2</sub> HPO <sub>4</sub>	6.99	0.025	+0.03
Borax	3.80 g Na <sub>2</sub> B <sub>4</sub> O <sub>7</sub> • 10 H <sub>2</sub> O	9.18	0.021	-0.05
Calcium hydroxide	5 g Ca (OH) <sub>2</sub>	12.48	0.13	

(\*) The pH value as a function of temperature is shown in the graph in fig. 3.1

NB. - Other solutions are given in the DIN 19266 paper and in IUPAC and IEC standards now being approved.

Fig. 3.1. pH values of the main standard solutions, as a function of temperature.

The degree of buffering  $B$  is defined as the inverse of the pH variation following the addition of an infinitesimal amount of a strong acid (HCL) or a strong base (NaOH):

$$B = \frac{1}{V_0} \frac{dn}{dpH}$$


where  $V_0$  is the volume of the solution (in litres),  $dn$  is the amount of added HCL or NaOH (expressed in equivalent grams) and  $dpH$  is the resulting pH variation.

The effect of dilution is defined as the pH variation resulting from the addition of water in a 1:1 ratio. The pH value of buffer solutions may be calculated by means of the mass action law, if we know the dissociation constant of the acid or the base, its concentration and the concentration of its salt. The dissociation constant may in turn be calculated by measuring its conductivity; since this depends on temperature, the pH of buffer solutions is also a function of temperature (**fig. 3.1**).

A certain amount of care is required in preparing buffer solutions. The substances used must be «analytically pure» and must be dissolved in demineralized or twicedistilled water with conductivity < 2 micro Siemens/cm. For the preparation of alkaline buffer solutions, water must be free from carbon dioxide, so it should be boiled before use. The duration of buffer solutions is not unlimited; if they are left in contact with the air for long periods, they deteriorate in that they tend to absorb CO<sub>2</sub> and there is a formation of mould.

## 4. REDOX POTENTIAL (ORP)

If electrons are given off during a chemical reaction this is known as oxidation; if on the other hand electrons are absorbed, this is known as reduction. It therefore follows that the substances that can accept electrons are oxidisers. A redox reaction (oxidation-reduction) may be written in the following general form:

When a noble metal is immersed in the reagent, a potential is formed whose value depends on the ratio between the concentration of the oxidiser and the concentration of the reducer. The more oxidising the reaction environment, the more positive the potential. Platinum is the most frequently used noble metal. Before a platinum electrode can be used it must be cleaned mechanically using abrasive substances, then immersed in diluted hydrochloric acid and finally washed with water. In some cases, e.g. if platinum has a catalysing effect or if the speed of response is slow (this also depends on the reversibility of the reaction), it will be preferable to use a gold electrode. Gold is not recommended in the presence of cyanide or chloride. The formation of air bubbles or other gases on noble metals should be avoided.

### 4.1 SPECIFIC ION ELECTRODES

Selective ion electrodes possess the following characteristics:

- 1) the signal is a voltage that must be measured without charging the electrode (high impedance).
- 2) the signal is proportional to the logarithm of the ionic activity.
- 3) the signal does not depend on the sensitive area of the electrode.

Ideally, a specific ion electrode should be sensitive only to one ionic species (similarly to the almost ideal situation of the pH electrode). In actual fact there is a more or less marked dependence on other ions (or other gases if the electrode measures a dissolved gas, for example ammonia). The degree of interference is specified by the manufacturers of selective ion electrodes. In many cases the solution must be titrated beforehand to an optimal pH value and/or ionic concentration. The solutions for preparing the specimens are called TISAB (total ionic strength adjustment buffer) or Condisol (conditioning solution). The measuring technique (preparation of specimens) influences the measurement result and therefore it must be specified. At 25°C a monovalent selective ion electrode (for example Na<sup>+</sup>) gives 59mV for every decade of concentration; for bivalent ions (e.g. Ca<sup>++</sup>), values of 29.5 mV decade are obtained. This means that an error in measurement of  $\pm 1$  mV leads to an error in terms of concentration of  $\pm 4\%$  with monovalent ions and  $\pm 8\%$  with bivalent ions. From the instrumental point of view it is possible to obtain greater precision than  $\pm 1$  mV, however this is not accompanied by a significant improvement of the final result, since the greatest limit is caused by the diffusion potential of the reference electrode. This potential is created at the ends of the porous septum and depends on ionic concentration. With special reference electrodes and/or by titration, uncertainty may be less than  $\pm 1$  mV.

## 5. CALIBRATION

The following adjustments are possible with the HD 8602:

- 1) compensation of electrode temperature
- 2) correction of electrode asymmetry ( $\Delta$  pH)
- 3) correction of electrode  $\sphericalangle$  slope.

**The correction of asymmetry potential must be carried out before adjusting the slope. Before calibrating the slope, the correct temperature compensation must be checked as described in paragraph 5.1.**

Temperature has relatively little influence on the correction of asymmetry potential (ideally it has no influence at all if the isothermal point coincides with pH = 7), so that it is possible to omit temperature compensation, unless a very high degree of precision is required (< 0.05 pH). For calibration we use solutions whose pH value is already known, called buffer solutions.

The calibration of the asymmetry potential (6.86 pH) should be checked every 10-20 days, according to the degree of precision required. Slope calibration (from pH 4.01 to pH 9.18) normally need only be carried out when the electrode has to be replaced.

## 5.1 COMPENSATION OF pH ELECTRODE TEMPERATURE

As mentioned in section 2, the voltage generated by a pH electrode increases linearly with the increase in temperature with a coefficient of  $\alpha = 0.335\%/^{\circ}\text{C}$ . To obtain the correct pH value, this dependence on temperature must therefore be compensated. With the HD 8602 it is possible to carry out automatic compensation and manual compensation. **Manual compensation is operative only if the temperature probe TP870 is not connected to the corresponding input on the instrument.** The procedure for manual compensation is the following:

- 1) Press button  $^{\circ}\text{C}$  to obtain the temperature reading.
- 2) Press button CAL to proceed to the calibration mode. The  $^{\circ}\text{C}$  symbol should flash.
- 3) Press button  $\Delta$  (to increase) or  $\nabla$  (to decrease) until the temperature shown is the same as that of the solution being analysed.
- 4) Press button CAL again to return to the measuring mode. In the measuring mode the  $^{\circ}\text{C}$  symbol does not flash.

With the temperature probe TP870 compensation is automatic, it is sufficient to insert the connector in the input provided on the instrument and immerse the probe in the solution in which the measurement is being taken.

## 5.2 CORRECTION OF ASYMMETRY POTENTIAL ( $\Delta$ pH)

Ideally a pH electrode supplies a signal of 0 V for pH = 7. In fact the electrode generates a potential that is not equal to 0 V (asymmetry potential) and therefore this potential must be compensated. The procedure is as follows:

- 1) Immerse the pH electrode (or the glass electrode and the reference one if you are using a combined electrode) in a buffer solution with pH = 6.86. If the probe for TP870 is available, immerse this too in the same solution. If this probe is not available, the temperature of the buffer solution must be set manually (par. 5.1). Switch the instrument ON, set the pH measurement.
- 2) Press button CAL to proceed to the calibration mode. The unit of measurement symbol on the display will flash. The symbol  $\Delta$  (for  $\Delta$  pH) will appear on the display. The instrument displays the pH value of the buffer solutions at the temperature measured by the probe, or manually set (automatic calibration).
- 3) Press button  $\Delta$  or  $\nabla$  until the display indicates the pH value of the buffer solution at the temperature taken. For example, set 6.86 pH if the buffer solution is at a temperature of  $25^{\circ}\text{C}$ , set 6.83 pH if it is at  $50^{\circ}\text{C}$ .
- 4) Press button CAL again to confirm the indicated pH value. The symbol  $\Delta$  disappears and the pH symbol stops flashing. If a very high degree of precision is not required ( $< 0.05$  pH), temperature compensation may be avoided during  $\Delta$  pH calibration (point 1). In fact the isothermal point,

that is the pH value for which the potential generated by the pH electrode does not vary as temperature varies, is close to pH = 7.

## 5.3 CORRECTION OF SLOPE

According to Nernst's law, an electrode supplies a variation of potential of 59.16 mV/pH at 25°C. A real electrode supplies a signal that differs slightly from this theoretical value, for which reason pH-meters must offer the possibility of adjusting the slope of the electrode.

The procedure is the following:

- 1) Immerse the pH electrode (or the glass electrode and the reference one if you are not using a combined electrode) in a buffer solution with pH = 9.18. If the probe TP870 for measuring temperature is available, immerse this too in the same solution. If this probe is not available, the temperature of the buffer solution must be set manually (cfr. par. 5.1).
- 2) Press button  and the button CAL. The pH symbol flashes.
- 3) Press buttons  $\Delta$  or  $\nabla$  until the display indicates the pH value of the buffer solution at the temperature taken. For example, set 9.18 pH if the buffer solution is at a temperature of 25°C, set 9.01 if it is at 50°C, etc.
- 4) Press button CAL again to return to the measuring mode. The symbol mV will appear on the display which will indicate the mV generated by the pH electrode. Press the pH button to obtain a reading of the pH value.

The calibration is completed.

## 6. MEASUREMENT

### 6.1 Measurement of the pH value

Connect the instrument to a pH electrode suitable for measuring requirements (e.g., use a pointed electrode for measurements in cheese and meat, a flat electrode for measurements on surfaces, etc.). Connect the temperature probe TP 870 and ensure that it is at the same temperature as the pH electrode. If the TP 870 probe is not available, set the temperature of the pH electrode manually. Note the pH value once the reading is stabilized.

If you stir the pH electrode about in the solution to be analysed you will observe a slight change in the reading. This phenomenon is due to a change in the diffusion speed of the ions which migrate from the solution inside the reference electrode into the solution to be analysed. If calibration is carried out while the buffer solution is being stirred (e.g. with a magnetic stirrer), then measurements must also be taken while the unknown solution is being stirred. The same applies to a still electrode. When the solution is stirred there is a quicker response and the readings are more easily reproduced. **When passing from one solution to another with a notably different pH value, there is the risk that residue of the first solution will cling to the electrode, thus threatening to alter the pH value of the**

**second solution.** Solutions with a low ionic concentration (low conductivity) are particularly prone to pollution. However we do not advise cleaning the electrode with paper towels, rags, etc. When the glass of the electrode is rubbed this generates an electrostatic charge that has a transitory worsening effect on reproducibility. The most practical way of cleaning the electrode is to divide the solution you wish to analyse into two parts: use the first part only for rinsing the electrode and take the measurement in the second part.

## 6.2 MEASUREMENT OF THE REDOX POTENTIAL AND MEASUREMENTS WITH SPECIFIC ION ELECTRODES

Connect the electrode for measuring the redox potential (ORP) or the specific ion electrode (ISE) to the BNC connector on the instrument. Press the mV button and read the value shown on the display. Specific ion electrodes are supplied with a chart or a graph showing the ionic concentration as a function of the mV generated. If the electrodes are not combined - as is frequently the case with specific ion electrodes - a reference electrode must also be connected to the input provided on the instrument with a Ø4 mm banana connector. During mV measurements, temperature does not influence the reading. In fact the temperature coefficient of specific ion electrodes depends on the valency of the ion, therefore it is not single. When measuring the redox potential it is impossible to compensate the temperature coefficient on account of the many interfering factors.

## 6.3 TEMPERATURE MEASUREMENT

Connect the temperature probe TP870 to the corresponding input on the instrument. Press the °C button and note the value given once the reading is stabilized. The TP870 probe is particularly suitable for immersion measurements. It is possible to measure the surface temperature of solids by spreading a paste with a high degree of thermal conductivity on the surface, for example silicone grease, and immersing the tip of the probe in the paste. Measuring the temperature of a gas is possible only as far as regards the average value over a period of a few minutes; this is due to the time constant which, in still air, can reach a value from 100 to 300 times higher than in a liquid.

## 6.4 MEASUREMENT OF THE pH OF EARTH

There are various methods of measuring the pH of earth. As the different methods give different results, it is essential to quote the method used as well as the pH value obtained. Table 6.4.1 shows the pH values obtained on immersing a combined pH electrode in a mixture of earth and distilled water.

**Table 6.4.1** - Effect of the earth/water ratio on the pH value

ml of water added to 100 g of earth	pH of the suspension		
	earth 1	earth 2	earth 3
10	7.45	9.10	7.95
25	7.60	9.40	8.00
100	7.70	9.85	8.20
1000	8.15	9.90	9.20

The glass electrode is immersed in the sediment (obtained by letting the suspension rest for a time), while the porous septum (or the reference electrode if this is physically separate from the glass electrode) remains in the liquid on top. More reproducible recordings may be made by measuring the pH of a suspension of earth in a neutral salt solution. Solutions of KCl 1M or CaCl<sub>2</sub> 0.01 M are generally used. These give pH values lower than those obtained for the same earth suspended in water;

**Table 6.4.2** - pH values measured in different conditions

earth suspended in:	pH recorded in:		
	liquid above	suspension	sediment
distilled water	6.2	5.8	4.7
KCl 1M	5.1	5.1	5.1

furthermore, the value measured is the same in the solution, the suspension and the sediment (contrary to what occurs when using a mixture with distilled water). For some types of earth, the earth-water mixture has a very low electric conductivity; in these cases the measurement is unstable and is not easily repeated. If on the other hand a solution of KCl is used, excellent repeatability may be obtained even in these difficult cases.

## 7. ERRORS IN MEASURING THE pH VALUE

There are basically two types of measurement error: instrumental errors and errors due to the sensor (electrode). Instrumental errors ( $\pm 0.05$  pH) are generally negligible in comparison with those due to the electrode. Especially for values of pH  $> 12$  and in the presence of  $\text{Na}^+$  ions there are deviations from the ideal behaviour due to the composition of the pH sensitive glass (about 0.1 pH). Another source of error is the diffusion potential created at the ends of the porous septum. For diluted solutions of acids, bases and salts with a concentration of less than 0.1 mole/litre there are errors lower than  $\pm 0.05$  pH; in a base or a concentrated acid there may be errors of  $\pm 0.2$  pH. These errors may be minimized by calibrating with solutions that have a pH value and an ionic concentration similar to that of the specimen. In specimens with suspended particles there may be errors of up to  $\pm 2$  pH. These errors may be eliminated by filtering the solution. In solutions with a very low electric conductivity, pH measurements are unstable and are not easily repeated; when a special solution is added (pHiX adjustor), the ionic concentration and consequently the conductivity are increased and the measurements become more stable. The change of pH due to the addition of the pHiX adjustor is negligible ( $< 0.003$  pH).

The glass membrane of a pH electrode must always be well hydrated.

If the electrode is immersed in an anhydrous, water-greedy substance (such as alcohol), its hydrated layer may be destroyed. Likewise this layer may be destroyed if the electrode is allowed to dry in the air. In contact with water the hydrated layer regenerates slowly, so the measurements are not reliable in the first two hours (passing from dry to damp).

**Fluorhydric acid, some fluorine compounds with pH  $< 5$ , strongly alkaline solutions at a high temperature and saline solutions in high concentration and at a high temperature may attack the glass membrane, so it is not advised to carry out measurements in these conditions.**

## 8. BUTTON FUNCTIONS

turns the instrument on or off. The instrument turns itself off automatically about 8 minutes after the last key was pressed. Automatic cut-off may be avoided pressing simultaneously HOLD and ON/OFF when switching on.

when this button is pressed the word HOLD appears on the display and the reading is stopped. When it is pressed for the second time the word HOLD disappears and the display is updated again.

when this button is pressed the instrument is set in the calibration mode and the buttons  $\Delta$  and  $\nabla$  are activated. The second function of the FUNCTION buttons (blue part) is also activated. The unit symbols flash on the display to indi-

cate that the instrument is in the calibration mode. When CAL is pushed for the second time the instrument returns to the measuring mode and the units stop flashing.


the function of these is to increase or decrease the value shown on the display if the instrument is in the calibration mode (the unit must be flashing).

for setting pH measurement with the instrument in the measuring mode; if the instrument is in the calibration mode, this activates calibration of the asymmetry potential (  ).

for selecting mV measurement with the instrument in the measuring mode; if the instrument is in the calibration mode, this activates  slope calibration.

if the temperature probe TP870 is connected to the instrument, the temperature of the probe is displayed. If it is not connected, the display shows the temperature set for the manual temperature compensation. When the key is pressed the display alternately shows the value in °C or °F.

## 9. CHANGING THE BATTERY

Press your thumb on the ridged part of the battery compartment cover and push towards the bottom of the instrument. Replace the battery and slide the cover back into the groove. When the flat battery symbol appears on the display there is still an autonomy of 30 minutes. After changing the battery both the instrument and the electrode have to be recalibrated.

## 10. MAINTENANCE

Regular maintenance of the pH electrode ensures optimum performance, longer life and less costs for part replacement. The instrument and the temperature probe require no maintenance. The pH electrode should be checked every 10-20 days for cracks, formations of salt crystals or deposits around the porous septum. Wash the saline deposits with distilled water. Various deposits around the porous septum and on the pH-sensitive membrane may be removed by immersing the electrode in a 0.1 molar solution of HCl or HNO<sub>3</sub> for 15 minutes.

A blocked porous septum may be cleared by immersion in KCl 0.1 M at 70°C (60°C for calomel electrodes) for 15 minutes. Empty the reference electrode, if possible, and fill it with fresh electrolyte. Immerse the electro-

de in a solution of KCl 0.1 M for 24 hours before use. If the porous septum remains blocked you must identify the substance causing the blockage and use a specific solvent. For example, protein deposits may be removed with 1% pepsin in 0.1 mole of HCl; grease and oil may be dissolved in a detergent or a solution of methyl alcohol. For an electrode to be always ready for use, it must always be well hydrated. We advise keeping the electrode in a storage solution or in a saturated solution of KCl. If no KCl is available, use a buffer solution with pH = 7. Otherwise the electrode may be preserved by covering it with a protective cap containing a saturated solution of KCl. Do not use demineralized water! The level of the solution that fills the reference electrode must be higher than the level of the solution in which the electrode is immersed (a few cm. higher). This avoids a flow towards the internal solution in the reference electrode. The electrode connector and the BNC input must not be wet. A wet or dirty connector may be cleaned by washing it in a dehydrating solvent that does not leave residue (e.g. alcohol). If the electrode is disconnected from the instrument, close the BNC input with the plug provided.

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### CONDITIONS OF GUARANTEE

All our electrical appliances have been submitted to accurate testings and are guaranteed for 24 months from date of purchase. The Firm will replace or at our option repair that part without any charge found defective within the guarantee period. This excludes the complete replacement of the appliance, provided that the appliance has been used in accordance with the instructions. All repairs will be carried out at our Technical Assistance Centre. The transporting charges are at the customer's expenses. The guarantee excludes: accidental breakage due to transporting, misuse, connection to an unsuitable electricity supply, probes, sensors, electrodes and all accessories. The guarantee also excludes damages due to neglect or repairing by personnel not authorized by us. The guarantee is operative only if the guarantee coupon has been correctly filled in all details and sent back to the Firm within 15 days from purchase date, the guarantee certificate retained by the purchaser must be shown to the authorized service agent or attached to the appliance should it need to be sent to one of our service centres. Any controversy falls within the cognizance of the Padova Judiciary Court.

**HD 8602**

Instrument type

Serial number

This equipment has been fully inspected, and conforms to the current product specification.

Inspector

Date

## RENEWALS:

Date

Date

Date

Inspector

Inspector

Inspector

## Return Information

Instrument type

**HD 8602**

Serial Number \_\_\_\_\_

Customer Address

Nature of Fault

Date

---

---

---

DELTA OHM SRL - VIA G. MARCONI, 5  
35030 CASELLE DI SELVAZZANO (PD) - ITALY

DELTA OHM SRL  
VIA G. MARCONI, 5 - 35030 CASELLE DI SELVAZZANO (PD) - ITALY  
TEL. 0039-0498977150 r.a. - FAX 0039-049635596  
e-mail: [deltaohm@tin.it](mailto:deltaohm@tin.it) - Web Site: [www.deltaohm.com](http://www.deltaohm.com)

NM HD8602

SAFIGRAF 03/02